

Church Puts the People with Hearing Loss in the Loop!

Proud supporter of



About 134,958 people in Jefferson, Oldham and Shelby Counties suffer from hearing loss. The church has installed a Hearing Loop System that will make hearing our service much easier for those with hearing loss.

Hearing loop systems take sound straight from the speaker's microphone and deliver it right into the listener's hearing aid without extraneous noise or blurring. To them, it

sounds like the speaker is right in their head. It turns their hearing aids into wireless earphones that broadcast sound customized for their hearing loss.

Hearing Loops are the only assistive listening system to send clear, pure sound directly to hearing aids. They are the international standard for universal hearing access.

How do you use the hearing loop?

All the user needs to do is switch his or her hearing aid to the "T-coil", telephone, or hearing loop program. All but the very smallest modern hearing aids tend to be supplied with an internal t-coil. However, it may not have been enabled when the device was initially programmed. If the user is not sure whether his or her hearing aid has a t-coil and/or does not know how to access the program, his or her audiologist or hearing aid specialist will be able to advise. Once your t-coil is set, you can use it anywhere that has a hearing loop.

Can a hearing loop be used by someone without a hearing aid?

Yes! We have several portable headsets available. They are not as convenient to use as a hearing aid with a built-in t-coil and are not fine-tuned to suit individual patterns of hearing, but can still be helpful to people with hearing loss who do not wish to wear hearing aids on a regular basis. Each receiver has a headset connected to a small control unit with a volume control knob. If you would like to use one of these units, ask one of the greeters to get you one. If you want to test the system out with your hearing aid or try the headsets during the week, call the church (502-245-0573) during business hours and make an appointment to get "in the loop".



Hearing Loop
Switch hearing aid to T-coil